

Patent Application of

Cleon L. Griffiths

for

TITLE: SECURE CONDOM

BACKGROUND-FIELD OF INVENTION

This invention relates to Condoms, specifically to an improved Condom that is designed to remain securely on male member until forcibly removed.

BACKGROUND-DESCRIPTION OF PRIOR ART

Condoms are made of very thin rubber or plastic in the form of tube that is closed on one end and open on the other. The tube is uniform in diameter for the entire length. The purpose of a Condom is to prevent unwanted pregnancies and prevent the spread of venereal diseases. The Condoms of prior art have several drawbacks. One drawback is the fact that if Condoms are used without lubrication inside, the sensation of intercourse is like going wading in a puddle with overshoes on. The advantage of using a condom dry is that it will likely remain on the penis thruout sexual intercourse. If a condom has lubrication inside, The sensation will be much better, but the condom is likely to come off during intercourse which would cancel out all of the benefits that you were using the condom for in the first place. ~~The~~ Secure Condom addresses both the problem of sensation with lubrication and an enlarged closed end for comfort, then it has the remainder of the condom too small to easily fit over the head of the particular penis it fits. (these come in various

COPIES OF THIS DOCUMENT ARE AVAILABLE FROM THE NATIONAL ARCHIVES AT COLLEGE PARK, MARYLAND

sizes and are packaged in lubricant.) Secure Condoms come with the small diameter area rolled and stretched upon an Installation ring. (fig 2) This ring is made of plastic or some other durable material. The outside circumference of the installation ring has a groove to accept the rolled condom. One side of the installation ring is of larger circumference to prevent the condom from being rolled off the wrong side of the ring. The inside of the installation ring is smooth and of large enough diameter to allow a penis to enter thru the ring into the enlarged area. (fig 3) To use this condom, it is removed from the package. The Secure condom is attached to the installation ring prior to packaging. Identifying the side of the ring that has the largest circumference tells the user that this is the front of the condom. The penis is pushed thru the installation ring into the enlarged area of the condom at this time the area of the condom rolled up on the ring is unrolled to the base of the penis. (fig 4) The installation ring is then slid over the head of the penis and discarded. (fig. 5) The area B of the condom is smaller in diameter than the penis that it fits, and will not go over the head of the penis without some effort. Sexual intercourse will not cause area B of the condom to go over the head of the penis. The fit on the shaft of the penis will be tight, but not uncomfortable. It is for this reason that the name of these condoms is SECURE CONDOMS BECAUSE THEY WILL NOT COME OFF DURING INTERCOURSE.. The second advantage to the Secure Condom is that the sensitive area of the Penis known as the Glans Penis is in an area of the condom that is not the least bit restrictive and is bathed in lubricant similar to sexual intercourse without a condom.

The Glans Penis can slide around in the enlarged area of the Secure Condom with the sensitivity of having sexual intercourse without a condom, and without the worry of the condom coming off or splitting as some too tight condoms are prone to do.

The following prior art patents have been found to be revelant to the field of the present invention.

US Pat. # 4,821,742 Issued to Pnelps on April 18, 1989 for "Contraceptive Device" Hereafter the "Phelps Patent"

US Pat. #4,966,166 issued to Leffler on Oct. 30 1990 for Prophylactic device and Method" Hereafter the "Leffler Patent"

US Pat. # 5,010,871 Issued to Christina on April 30, 1991 for "Prophylactic Safety Device" Hereafter the "Christina Patent"

US Patent No. 5,111,831 issued to Foggia on May 12, 1992 for "Scrotum Supporting Condom with Retention Means" Hereafter the "Foggia Patent".

US Patent No. 5,199,444 issued to Wheeler on April 6, 1993 for " Condom having enhanced Grippability structure and annular sealing element" Hereafter the "Wheeler Patent".

US Patent No. 5,284,159 issued to Wilk on Feb 8, 1994 for "Prophylactic Device" Hereafter the "Wilk Patent".

The Phelps Patent is a cap that is attached to the head of a penis with adhesive. The main drawback is the pain of removing this cap from a very sensitive area. Also the penis is very seldom completely dry and would not hold adhesive very well so the possibility of the cap coming off during intercourse is great. This is one of the better features of the Secure Condom is that it will stay on a penis during intercourse without the use of adhesive or straps or other contrivances.

The Leffler Invention has the advantage of having more sensitivity than regular condoms and would probably stay on, but

Page 4

it is complicated to make and use. There is a statement in the claims that shows a definite weakness in the use of this condom. "The device of claim 1, wherein the inside of the mid-section of the device is lubricated and wherein the inside edge of the sleeve section is closed and releasably sealed to itself thereby separating the lubricated inside of the midsection from the inside surface of the sleeve section." This problem of a condom being lubricated in part but not everywhere would make it impossible to lubricate it before packaging without the lubricate getting all over the condom. There could also be a problem with installing this condom thru these sleeves. With the Secure Condom the manufacture of it would be no more difficult than a regular condom of prior art. The Secure Condom can be packaged in water soluble lubricant so that it is lubricated all over and not just in a certain area, which makes it much simpler to make. With the smaller circumference of the Secure Condom that fits the shaft of the penis tightly the user will have the advantage of having the sensitivity of lubrication inside and outside with no worry of the condom coming off or leaking.

US Patent No. 5,010,871 issued to Christina on April 30, 1991 for "Prophylactic Safety Device" This is a ring that is apparently slid on the condom and onto the base of the penis. It is packaged with this ring attached to a corner of the open end of the condom by a strip of similar material as the condom. This Condom does not address the movement of the head of the penis in an enlarged lubricated area such as the Secure Condom does, but only keeps the condom from coming off. Also this ring could be uncomfortable and cut off blood flow to the penis if too tight. Also all men do not have the same size penis, so even though this

ring was comfortable for one user, a user with a larger penis would find it uncomfortable. It is also more complicated in use than the Secure Condom.

Foggia Patent; This patent of a condom that has some rubber parts that encircle the scrotum as well as an inner seal of smaller diameter than the rest of the condom. Otherwise the configuration is the same as regular condoms. This condom is difficult to manufacture due to the contrivances attached to a regular condom. This condom could be difficult to use in respect of pushing the penis thru the small seal, and very uncomfortable having a strap around the scrotum, even one of the same material as the rest of the condom. The Secure Condom does everything that this invention does and more without some awkward contrivances to contend with. It will also be much cheaper and easier to produce.

Wheeler Patent; This patent has a flap of similar material to the condom glued or otherwise attached at the open end of the condom with a smaller opening than the remainder of the condom which otherwise resembles a regular condom of prior art. In this invention one embodiment has the rest of the condom rolled up on an applicator ring, however the smaller circumference compressive retention means is not stretched. The drawing of this, fig. 21 and 22 on his patent shows this smaller area not stretched onto the ring. He doesn't say how he got the rest of the condom on the ring without getting this part on, but that must be a mystery. Anyway it would be difficult to install this condom by forcing this small opening over the penis of the user. With the Secure condom the smaller diameter area of the condom is completely rolled upon the installation ring so that the penis can enter the enlarged area

at the closed end of the condom, making it very easy and effortless to install. With the entire area of the shaft of the penis being tightly fit with the Secure condom, not only will this prevent leakage, but will keep the condom on more securely. These condoms will be packaged in soluble lubricant which will make the production of the Secure Condom much simpler to make but also simpler to use. One embodiment of Wheelers invention is having a strip at the open end that is covered with adhesive so that the condom can be overlapped to make it tight at the open end. Again this condom would be difficult if not impossible to have lubricant packaged along with the condom or the lubricant would get on the adhesive strip off the hands when trying to attach the adhesive to itself thereby causing the adhesive to not stick cancelling this benefit. If the condom is not lubricated, you have lost the benefit of the sensitivity of having lubrication in the condoms. The Secure Condom has none of these limitations. It is much simpler to produce than the above inventions. It can be lubricated inside and outside when it is packaged. It can be installed very simply even in the dark since the installation ring will have a larger circumference on the front side which will make it impossible to roll it off the wrong way. It has an enlarged area which will give the head of the penis the sensitivity of not even having a condom installed and it will not come off by accident if the correct size Secure Condom has been chosen for the particular size penis that it is used on. These condoms will come in different sizes to fit different size penis. The main advantages of Secure Condoms is their simplicity to produce, and to use. Another main advantage is the fact that

lubricated inside and outside when packaged. Another big advantage is the enlarged area at the closed end of the Secure Condom that lets the sensitive head of the penis slide freely within the condom during intercourse without fear of it coming off. The sensitivity enjoyed by the user of Secure Condoms is as near as it is possible to achieve when using Condoms.

There is a desperate need for condoms to prevent disease and unwanted pregnancies, however the prior art leaves much to be desired. When the condoms are an improvement over old type condoms, yet the methods that they have come up with are not economical to produce, not practical in use or too difficult to use. The Secure Condoms design is economical to produce because they will be manufactured using current procedures. The Secure Condom will be welcomed because it is lubricated inside and out to give as much sensitivity as is possible in a condom. The enlarged area will give a feeling of freedom that current condoms do not have. In fact Secure Condoms will feel to the user as though he does not have a condom on at all. The Secure Condoms will have the smaller diameter section of the Condom rolled upon the installation ring which is grooved to hold the condom. This installation ring will have one side that is of larger diameter than the other side of the groove to indicate that this side is the front of the condom, thus it will be easy to install even in the dark. The enlarged area of The Secure Condom will not be rolled up on the installation ring but will extend thru the ring. There will be no interference of the part of the condom rolled up on the ring so that a penis can freely enter the enlarged area thru the installation ring. To use these Secure Condoms. 1. Remove from package. 2. find

Page 8

which side of the ring is larger; this is the front. 3. push your penis thru the ring into the enlarged area. 4. Roll the part of the condom that is on the ring, off the ring and all the way back to the base of the penis. 5. slide the installation ring off over the nose of the penis and discard.

The Secure Condom is now installed. Wilks patent is the same as contemporary state of the art condoms. The only difference is that the Wilks Patent is double thick and has holes at the open end to equalize the pressure of the outer and inner condoms. There is nothing about this Condom to aid in preventing the condom from coming off during intercourse if it is lubricated on the inside for better more natural sensation. The Secure Condom has several advantages over the Wilks Patent. The Secure Condom has an enlarged area for the first 2 inches at the closed end which allows movement inside the condom of the Glans Penis, which is the sensitive area of the Penis. The Secure Condom will be packaged in water soluble lubricant so it will be lubricated inside and out for a more natural sensation. Since it will have a smaller area behind the enlarged area that will not go easily over the head of the penis, there is no fear that the secure condom will come off during intercourse. The Secure condom will have a ring upon which the smaller circumference area of the condom will be rolled. This ring will be made of plastic or other strong material and will have an inside diameter large enough for the penis will be able to enter the enlarged area of the condom unimpeded. After the penis is pushed thru the ring, the part of the condom that is rolled up on the ring will be rolled off the ring and unrolled all of the way to the base of the penis. At this time the ring will be slid off over the head of the penis and discarded. The Secure Condom is now installed.

Page 9

Secure Condoms will come in different sizes, since all men are not endowed equally. Cardboard gauges will be available at no charge, so a man can determine what size condom is right for him. The ring for each size condom will be the same inside diameter as one of the gauges. The condoms will also come in different sizes to match the gauges.

SECURE CONDOMS
SUMMARY

The Secure Condom addresses all of the weaknesses of condoms of prior art. A list of these weaknesses follows.

1. condoms accidentally coming off during sexual intercourse.
2. condoms splitting from being too tight.
3. condoms coming off if lubricant is used.
4. lack of sensitivity due to being too tight around the sensitive area known as the Glans Penis.

Secure condoms come in various sizes to fit male members of various sizes. The Secure Condom fits the male members quite tight from the base of the penis to the rear of the head of the penis so that some effort is required to remove this condom. over the head of the penis, therefore it will not come off due to the activity of sexual intercourse. (1), (3)

Secure Condoms have an enlarged area for about two inches at the closed end of the condom which is somewhat larger than the head of the penis. This enlarged area gives a sensation of not having a condom installed, and in conjunction with the fact that these condoms are lubricated inside and out, the sensation is the same as no condom at all. (2), (4)

Secure condoms will be packaged rolled up onto the installation ring with water soluble lubricant. Since there is no fear of these condoms coming off, being lubricated inside and out is no problem, and with the enlarged area at the closed end there is no fear of splitting. (1), (2), (3), (4)

Secure condoms will come in various sizes since male members come in various sizes. Cardboard gauges with holes similar to

the holes in the different sized installation rings will make it easy to determine the correct size condom for each user. . . . these gauges will be free and available wherever secure condoms are sold. The installation rings will have the condoms small circumference area rolled up on the groove on the outside of the ring. The enlarged part of the condom will be pushed thru the ring before packaging. The installation ring will have the front side of the groove of larger circumference to prevent the condom from being rolled off in the wrong direction.

The manufacture of Secure Condoms will use the same methods of manufacture as the condoms of prior art most in use at present. Other than the cost of the installation ring, there is no reason that these condoms will be more expensive than the condoms of prior art.